

Ferrara Treatment Chair

"Our new Ocura Ferrara chairs can be raised to the perfect height for treatments which helps staff be more efficient and improves the patient experience."

Key Features

Our Ferrara chair is a compact treatment chair suitable for a wide range of treatments. It is available with upholstered cannulation type armrests or standard armrests, and can be easily placed in a comfort reclined position, or a flat "bed" position, making it a versatile therapy chair; it can also be put into an emergency "Trendelenburg" position quickly. With a seat width of 600mm and the perpendicular angle of the leg section, it can be used as a chair or couch. It is available with or without height adjustment and offers a wide range of accessories to make the Ferrara a next-generation chair.

- Easy to clean antimicrobial vinyl to aid infection control
- The vertical legrest enables the chair to be used in spaces that would be too small for a full length couch
- Innovative armrest technology with individual positioning
- Easy side access
- 100mm twin wheel castors
- (220kg) 34.6 stone

View a video of the Ferrara in action



BE1275
Ferrara treatment chair, in the Trendelenburg position



BE1275
Ferrara treatment chair, in the flat position, shown with accessories

Colour Options

Available in a wide range of colours, from our Valencia range, these can be found on page 206.





BE1270

Dimensions

All Models

Depth	1100mm
Seat Depth	500mm
Reclined Length	1960mm
Width	780mm
Seat Width	600mm

Height

Seat Height

Footrest dimensions Armrest size 1420mm
1420-1630mm (height-adjustable)
540mm
540-750mm (height-adjustable)

440mm x 130mm

550mm x 120mm (upholstered) 440mm x 50mm (standard)



Ordering Information

When ordering, quote the product code and add the colour code as a suffix, i.e. (Code)**AT** for Atlantic.

Code	Description
BE1260	Ferrara treatment chair, 3 motor, fixed height, upholstered arms, individually locking castors
BE1261	Ferrara treatment chair, 3 motor, fixed height, upholstered arms, central locking castors
BE1262	Ferrara treatment chair, 3 motor, fixed height, upholstered arms, fold-out footrest, individually locking castors
BE1263	Ferrara treatment chair, 3 motor, fixed height, upholstered arms, central locking castors, fold-out footrest
BE1264	Ferrara treatment chair, 4 motor, height-adjustable, upholstered arms, central locking castors
BE1265	Ferrara treatment chair, 4 motor, height-adjustable, upholstered arms, central locking castors, fold-out footrest
BE1270	Ferrara treatment chair, 3 motor, fixed height, standard arms, individually locking castors
BE1271	Ferrara treatment chair, 3 motor, fixed height, standard arms, central locking castors
BE1272	Ferrara treatment chair, 3 motor, fixed height, standard arms, fold-out footrest, individually locking castors
BE1273	Ferrara treatment chair, 3 motor, fixed height, standard arms, central locking castors, fold-out footrest
BE1274	Ferrara treatment chair, 4 motor, height-adjustable, standard arms, central locking castors
BE1275	Ferrara treatment chair, 4 motor, height-adjustable, standard arms, central locking castors, fold-out footrest

Accessories	
BE1280	Relax pillow
BE1281	Battery back-up
BE1282	Push handle with holder for paper rolls
BE1283	Foot switch for emergency CPR positioning
BE1284	One piece moulded polyurethane armrests (for upholstered arm versions only)
BE1286	Transparent cover for leg rest (not available on versions with fold-out footrest)
BE1287	Seat with pressure relieving springs
CA1825	Holder for infusion pole
CA1826	Infusion pole
CA1862	Oxygen bottle holder with assembling bridge

ocura.co.uk t: 01732 367777 e: sales@ocura.co.uk