



Mauro Relax Chair

“For one of our covid patients, being able to use one of these chairs made an enormous difference to his recovery progress and sense of wellbeing”



Dimensions

All Models

| | |
|-----------------|---------------|
| Width | 720mm |
| Depth | 840mm |
| Reclined Length | 1550mm |
| Seat Width | 535mm |
| Seat Depth | 510mm |
| Armrest Width | 75mm |

BE2062

| | |
|-------------|---------------|
| Height | 1290mm |
| Seat Height | 510mm |

BE2061/BE2064

| | |
|-------------|---------------|
| Height | 1210mm |
| Seat Height | 440mm |

BE2060/BE2065/BE2073

| | |
|-------------|---------------|
| Height | 1270mm |
| Seat Height | 510mm |


View a video of the Mauro Recliner in action



Key Features

The Mauro relax chair is equipped with a relax mechanism that tilts the backrest, seat and legrest into a zero-gravity position inspired by NASA*

The seat shell was designed following the shape of the human spine (S-shape) offering exceptional support.


- Easy to clean antimicrobial vinyl and a removable seat cushion to aid infection control
- Ergonomic design offering exceptional support for the lower back and the shoulders
- Efficient relax movement reducing pressure points on the body providing good blood circulation
- Large castors for easy movement (on selected models)
- The ergonomic armrests are height-adjustable and support 3 height positions, at their lowest position they are level with the seat allowing easy side access for patients with mobility issues
- Side handles can lock the recline into any position for safe use by the patient (applies to manual recline models)
-  23.5 stone/150kg
(Certain versions can be upgraded to 200kg, additional cost, please enquire)



BE2062 – 300mm front castors for patient transport







Colour Options

Our nine most popular colours are below, however these chairs are also available in the full Valencia range shown on page 206.




Ordering Information

When ordering, quote the product code and add the colour code as a suffix, i.e. (Code)AT for Atlantic.

| Code | Description | |
|--------|---|---|
| BE2060 | Mauro Ergo-Line – synchronised manual movement of backrest and legrest hospital patient reclining chair |  |
| BE2061 | Mauro Ergo-Line – synchronised manual movement of backrest and legrest hospital patient reclining chair with housekeeping castors |  |
| BE2062 | Mauro hospital transfer chair with 300mm front castors |  |
| BE2064 | Mauro Duo-Line – manual independently adjustable backrest and legrest hospital patient reclining chair with housekeeping castors |  |
| BE2065 | Mauro Duo-Line – manual independently adjustable backrest and legrest hospital patient reclining chair |  |
| BE2073 | Mauro electric backrest and legrest adjustment (independently or synchronised) hospital patient reclining chair |  |

All chairs have 4 x 125mm castors unless specified otherwise, *65mm housekeeping castors not suitable for patient transfer.

Models marked with  are available from our 48 Hour Express Range, in Atlantic.

| Accessories | |
|-------------|---|
| BE2013 | Mauro padded arm covers (pair) |
| BE2906 | Neckrest cushion |
| BE2021 | Mauro and Avento pressure care seat cushion |
| BE2063 | Anatomical reduction cushion |
| BE2067 | Mauro fully upholstered arm |
| BE2071 | Mauro/Avento detachable headrest |
| BE2076 | Extra wide swing out armrests |

*"Zero-Gravity" position – This sitting/lying position is inspired by NASA data on a neutral physical position, which transforms all muscular strength into a biomedical equilibrium



Mauro in Maternity

“The new Mauro reclining chair has allowed me to feel so much closer to my baby even though she has tubes attached, I can smell and cuddle her”

Mauro has improved the lives of parents and infants in maternity units, the goal of Mauro is to make kangaroo care, expressing and bonding accessible and comfortable for all recovering mothers and newborns. The recliner can be locked in the recline in any position, meaning the mum or partner feels safe holding their baby securely in a way that is comfortable.

The zero-gravity recline position means mum or partner can maintain kangaroo care, cuddles and feeding for an extended period without tiring. The armrests are height adjustable, at their lowest height they are flush with the seat cushion, allowing easy side transfer supporting new mums who have experienced a difficult birth or a C-section.

[Read more on how Mauro performs in maternity and neonatal](#)



