

Verona Therapy Chair

“I believe our new chairs have made a positive difference to the patients attending the thyroid biopsy clinic, and to the consultants examining the patients during this clinic”

Verona Therapy Chair – Key Features

Verona provides maximum clinical flexibility and the highest level of comfort in areas such as chemotherapy and dialysis. The Verona has a range of movement and the ability to go from a fully upright position to a fully reclined Trendelenburg position, allowing hospitals to use it for initial consultations and treatment. Equipped with up to five electric adjustable motors, the Verona is simplistic in its usability for both patients and clinicians with a picturised handset, to guide users into the desired position.

In addition, the polyurethane arms are cushioned and shaped to facilitate cannulation, and can fold back for easy access from the side of the chair. Complete with an optional shroud over the base and side of the chair, the Verona is easy to clean.

- Antimicrobial vinyl to aid infection control
- Transforms from upright to emergency Trendelenburg position in 7 seconds
- Independent motors allow for adjustments to the leg rest, footrest, seat tilt, backrest and seat height depending on model chosen
- Low entry height of 540mm ideal for patients with mobility issues
- Relax pillow supplied in light grey vinyl

-  34.6 stone/220kg





CA1815
Verona deluxe therapy chair in relax position



CA1815
Verona deluxe therapy chair in Trendelenburg position



CA1815
Verona deluxe therapy chair in bed position

Dimensions

All Models

Depth	1480mm	Height	1420mm
Seat Depth	500mm		1420-1630mm (height-adjustable)
Reclined Length	2300mm	Seat Height	540mm
Width	800mm		540-750mm (height-adjustable)
Seat Width	600mm	Footrest dimensions	2300mm x 800mm
		Armrest size	570mm x 130mm

Colour Options

Available in a wide range of colours, from our Valencia range, these can be found on page 206.



Ordering Information

When ordering, quote the product code and add the colour code as a suffix, i.e. (Code)**AT** for Atlantic.

Code	Description
CA1803	Verona therapy chair, 3 motor, fixed height, individually locking castors, manual footrest
CA1813	Verona therapy chair, 3 motor, fixed height, central locking castors, manual footrest
CA1804	Verona premium therapy chair, 4 motor, adjustable height, individually locking castors, manual footrest
CA1814	Verona premium therapy chair, 4 motor, adjustable height, central locking castors, manual footrest
CA1805	Verona deluxe therapy chair, 5 motor, adjustable height, individually locking castors, electric footrest
CA1815	Verona deluxe therapy chair, 5 motor, adjustable height, central locking castors, electric footrest

Accessories	
CA1820	LED reading lamp
CA1821	Protective cover – leg section
CA1822	Protective cover – footrest
CA1823	Battery back-up
CA1824	Push handle with holder for paper rolls
CA1825	Holder for infusion pole
CA1826	Infusion pole
CA1827	Foot switch for emergency Trendelenburg positioning
CA1828	Padded vinyl armrests
CA1862	Oxygen bottle holder with assembling bridge
BE1285	Premium Dolce C5 vinyl for Verona and Ferrara chair