




 New


curaPRO Vico Stretcher High/Low

Numerous possibilities and accessories make the work of nursing staff easier and improve patient welfare.

Key Features

- The practical handles help you manoeuvre the stretcher. When not in use, the push handles can be folded with ease. The push handles come as standard at the head end or they can be installed at the foot end
- Efficient to manoeuvre in small rooms and narrow corridors thanks to the 5th castor.
- The high/low adjustment is easy thanks to the hydraulic pump and foot pedal
- Easy-to-clean removable synthetic base cover aids infection control
- The backrest can be easily adjusted by to 90° by means of a gas spring
- The chair corners include a bumper to protect the chair and its surroundings from damage
- The side guards protect the patient over 2/3 of the stretcher
-  39 stone/250kg



Colour Options



Black

Ordering Information

Code	Description
BE2083	curaPRO Vico Stretcher High/Low

Accessories	
BE2084	2 Section foldable and adjustable IV rod at head end
BE2085	2 push handles at foot end
BE2086	Fixed oxygen bottle holder epoxy attached to the carriage
BE2087	Brake pedals at head end
BE2088	Mattress 10cm (with 3cm visco foam)
BE2089	Integral castors 200mm with central brake with 5th castor (antistatic directional castor)
BE2090	Mattress extension visco foam
BE2091	X-Ray Cassette Holder for Visco



Dimensions

All Models

Depth Upright	2175mm
Length in flat/bed position	2175mm
Width	800mm
Seat Width	690mm
Height	460-855mm