

**Flammability Rating:**

BS5852 Ignition Source (Crib) 5

**Antimicrobial Protection:**

Spectrum Colour Palette vinyl conforms to SNV 195 921:1994 for antimicrobial properties and is safe for use in the healthcare industry.

**Cleaning:**


Wipe clean routinely with a microfibre cloth, lightly dampened with cold water only, then dry with a clean, absorbent cloth. Remove body fluids and extreme soiling immediately (within 15 minutes) with cold water. Apply your approved infection protocol, a 5% sodium chlorite bleach solution, Haz-Tabs or Chlor-Clean. Remove product residues and then dry with a clean, absorbent cloth.

**Stains:**

Dyes and pigments from jeans, ballpoint and felt tip pens may permanently stain.

**Colours:**

Although every effort is made to ensure that the vinyl supplied matches the colours on this page, variations may occur due to limitations in printing processes.




Tin - TI

Grey - GY

Black - BC

Marine - MR

Blue - BL



Lagoon - LA

Sky - SK

Sea - SE

Grape - GR

Lilac - LI



LE MANS Electric Blood Sampling Chair in Grape

BORDEAUX ENT Ophthalmic Chair in Blue